

POSTER PRESENTATION

Open Access

Non-ECG-triggered perfusion imaging with integrated T1 mapping for quantifying myocardial blood flow

David Chen^{1,2*}, Behzad Sharif², Janet Wei³, Reza Arsanjani³, Xiaoming Bi⁴, C Noel Bairey Merz⁵, Daniel S Berman³, Debiao Li^{2,6}

From 17th Annual SCMR Scientific Sessions
New Orleans, LA, USA. 16-19 January 2014

Background

The dual bolus method is conventionally used for quantification of myocardial blood flow (MBF). The method necessitates an additional scan for each first-pass measurement, which increases total scan time and is sensitive to change in heart rates between the two scans. A single scan method to measure the arterial input function (AIF) using highly constrained back projection (HYPR) reconstruction for integrated T1 mapping was previously developed [1]. Mistriggering due to arrhythmias, changing heart rates, and poor electrocardiogram (ECG) signal may cause motion artifacts in HYPR reconstructed images because data is shared across multiple cardiac cycles. We propose using a non-ECG-triggered acquisition with cardiac motion self-triggering for integrated T1 mapping derived AIF for MBF quantification.

Methods

210 measurements were acquired continuously with no ECG triggering. 60 projections were acquired immediately following a saturation recovery (SR) magnetization preparation with a golden angle trajectory. Cardiac motion self-triggering images were produced from 60 projections using non-Cartesian SENSE. A sliding window was used to achieve a temporal resolution of 43 ms. The mean signal intensity in a region of interest drawn over the heart was measured. Peaks corresponding to increased ventricular blood pool volume in diastole were used for triggering. T1 maps were produced by measuring longitudinal magnetization recovery following a





Figure 1 Non-ECG triggered images of each TI during a single cardiac cycle. Images were created from 15 projections and used for T1 mapping.

¹Biomedical Engineering, Northwestern University, Evanston, Illinois, USA
Full list of author information is available at the end of the article



saturation recovery preparation [1]. Contrast agent concentration in the blood pool were found using T1 values. MBF was found using linear time invariant model-independent deconvolution. Six healthy volunteers underwent rest perfusion MRI studies on a Siemens 3T Verio. First, a dual bolus protocol was performed using a conventional clinical Cartesian sequence [2]. 10 min following the first perfusion scan, a second first pass perfusion experiment was performed using the proposed method. A two sided Student t-test was performed to MBF at a $p = 0.05$ significance level.

Results

Mean MBF found from the dual bolus and proposed non-ECG triggered, integrated T1 mapping acquisitions were 0.81 ± 0.29 ml/min/g and 0.93 ± 0.21 ml/min/g. There was no significant difference ($p = 0.43$) between MBF found using dual bolus and proposed acquisitions.

Conclusions

Non-ECG-triggered perfusion MR with fast T1 mapping for integrated AIF measurement produces similar MBF as gold standard dual bolus method. Fast T1 mapping alleviates the need for a separate scan specifically used for the acquisition of the AIF. Non-ECG-triggered imaging may improve robustness to arrhythmias and changing heart rates. This method may improve clinical feasibility of quantitative myocardial perfusion imaging.

Funding

NIH Grant T32 EB51705. NIH Grant RO1 EB002623. AHA Postdoctoral Fellowship Award 1POST7390063.

Authors' details

¹Biomedical Engineering, Northwestern University, Evanston, Illinois, USA. ²BIRI, Cedars-Sinai Medical Center, Los Angeles, California, USA. ³S. Mark

Taper Foundation Imaging Center, Cedars Sinai Medical Center, Los Angeles, California, USA. ⁴Siemens Medical Solutions, Los Angeles, California, USA. ⁵Women's Heart Center, Cedars-Sinai Heart Institute, Los Angeles, California, USA. ⁶Bioengineering, University of California, Los Angeles, Los Angeles, California, USA.

Published: 16 January 2014

References

1. Chen, Li: *MRM* 2013.
2. Ishida, Nagel: *JCMR* 2011, 13-28.

doi:10.1186/1532-429X-16-S1-P202

Cite this article as: Chen et al.: Non-ECG-triggered perfusion imaging with integrated T1 mapping for quantifying myocardial blood flow. *Journal of Cardiovascular Magnetic Resonance* 2014 **16**(Suppl 1):P202.

Submit your next manuscript to BioMed Central and take full advantage of:

- Convenient online submission
- Thorough peer review
- No space constraints or color figure charges
- Immediate publication on acceptance
- Inclusion in PubMed, CAS, Scopus and Google Scholar
- Research which is freely available for redistribution

Submit your manuscript at
www.biomedcentral.com/submit

